

Zebra ISBT 128 Labeling Solution



ICCBBA-Approved for ISBT 128 Labels

Blood is a precious resource. With Zebra's high-quality ISBT 128-compliant bar code blood bag labels, you can clearly label and accurately track every blood product—avoiding the waste and expense of having to destroy unidentifiable blood. You'll also reduce the potential for blood transfusion errors—a preventable error for which the Centers for Medicare & Medicaid Services will no longer reimburse.

One of just a few label converters licensed by the ICCBBA for ISBT 128 labels, Zebra delivers long-lasting, extremely legible blood bag labels. You're assured that once blood bags are labeled, they stay clearly labeled—including on blood products stored up to 10 years.



- Zebra labels are made of durable polypropylene material that is resistant to moisture—the label will remain intact whether you store it in a refrigerator or freezer.
- Designed to work together to provide optimal quality and durability, Zebra's thermal transfer blood bag labels and thermal bar code printers produce sharp, smear-proof, long-lasting bar codes, text and graphics that are resistant to moisture and common solvents—including alcohol and hand sanitizers.
- Zebra's blood bag labels adhere throughout the blood bag's lifespan. With their strong adhesion and pliability, they can conform to the shape of the blood bag to further ease scanning and ensure high bar code read rates.
- Our blood bag labels are available in a full-face configuration and are available in custom sizes.
- Zebra 8000T Primary Blood Bag label uses FDA 175.105-compliant adhesives, so it can be applied directly to the bag.

For more information about Zebra's blood bank labeling solution, call +1 800 423 0442.



Zebra ISBT 128 Label Media

- 8000T Primary Blood Bag label—Thermal transfer top-coated polypropylene label with a permanent acrylic adhesive that complies with FDA 175.105.
- 8000T Blood Bag label—Specialty thermal transfer polypropylene label with a permanent acrylic adhesive. Intended for use as an overlay only.
- Recommended genuine Zebra™ ribbons: 5555, 3200, 5586

Zebra® Easy-to-Use Printers

- ZM400™—Durable tabletop printer provides fast processing and printing speeds—making it ideal for collection centers that process a high volume of blood.
- G-Series™—Space-saving desktop printers are ideal for collection centers that process a low to medium volume of blood. Plus, they're lightweight and easy to carry on blood drives.



www.zebra.com

Corporate Headquarters
+1 800 423 0442
E-mail: inquiry4@zebra.com

Asia-Pacific Headquarters
+65 6858 0722
E-mail: apacchannelmarketing@zebra.com

EMEA Headquarters
+44 (0)1628 556000
E-mail: mseurope@zebra.com

Latin America Headquarters
+1 847 955 2283
E-mail: inquiry4@zebra.com

Other Locations

USA: California, Georgia, Rhode Island, Texas, Wisconsin **Europe:** France, Germany, Italy, Netherlands, Poland, Spain, Sweden **Asia Pacific:** Australia, China, India, Japan, South Korea **Latin America:** Argentina, Brazil, Florida (USA), Mexico **Africa/Middle East:** Russia, South Africa, United Arab Emirates

©2010 ZIH Corp. All product names and numbers are Zebra trademarks, and Zebra and the Zebra head graphic are registered trademarks of ZIH Corp. All rights reserved. All other trademarks are the property of their respective owners.

GSA#: GS-35F-0268N
P1032774 (11/10)

